



microman[®] **E**

Trade-in and get 30% off

Ergonomic pipetting for problem liquids

Specially created for scientists who work with problem liquids, the new MICROMAN E resolves common pipetting problems. With ease-of-use and control through an **ergonomic design**, a **patented QuickSnap** feature for Capillary Piston fitting and a **volume control button** to secure volume setting, MICROMAN E ensures accuracy, precision and comfort.

- Experience complete and rapid pipetting of viscous, dense, hot or cold samples.
- Avoid leaks when pipetting volatile (high vapour pressure) liquids.
- Avoid cross-contamination or residual sample carry-over.

Model	Volume µl	P/Number	Offer Price
M10E	1-10	FD10001	€278.50 €194.95
M25E	3-25	FD10002	
M50E	20-50	FD10003	
M100E	10-100	FD10004	
M250E	50-250	FD10005	
M1000E	100-1000	FD10006	

HOW TO TRADE-IN

Trade-in any pipette and get 30% off a new MICROMAN E. Email: suzanne@apexscientific.ie call **Dublin:** 01 6854686 or **Cork:** 021 2429776 for details. **QUOTE OFFER CODE: GIL095WE**

30% off Capillary Pistons

MICROMAN Capillary Pistons (CPs) offer 100% protection against contamination. Available as standard CPs or sterile CPs which are free of DNA, RNA, DNAses, RNAses, ATP, pyrogens and trace metals. Now with 30% off list price and new extended Sterile CP range.

HOW TO ORDER

QUOTE OFFER CODE: GIL094WE

- Call **Dublin:** 01 6854686 or **Cork:** 021 2429776
- Email: suzanne@apexscientific.ie
- Visit: www.apexscientific.ie/gilson

Capillary Piston	MICROMAN Model	Part Number	Qty / Box	€ Pricing	
				€ List	€ Offer
CP10	M10	F148412	192	€66.95	€46.87
		F148312	960	€266.60	€186.62
CP25	M25	F148112	200	€56.24	€39.37
		F148012	960	€293.33	€205.33
CP50	M50	F148113	200	€56.24	€39.37
		F148013	960	€293.33	€205.33
CP100	M100	F148414	192	€66.95	€46.87
		F148314	960	€266.60	€186.62
CP250	M250	F148114	200	€56.24	€39.37
		F148014	960	€293.33	€205.33
CP1000	M1000	F148560	182	€66.95	€46.87
Sterile Capillary Piston					
CP10ST	M10	F148413	192	€76.34	€53.44
		F148313	960	€378.62	€265.03
new CP25ST	M25	F148712	576	€227.17	€159.02
new CP50ST	M50	F148713	576	€227.17	€159.02
CP100ST	M100	F148415	192	€76.34	€53.44
		F148315	960	€378.62	€265.03
new CP250ST	M250	F148714	576	€227.17	€159.02
CP1000ST	M1000	F148180	182	€76.34	€53.44

Offers valid until 30 April 2015 to Ireland Customers only. *Terms & conditions apply



APEXScientific Limited | Unit F12 | Maynooth Business Campus |
Maynooth | Co. Kildare | Ireland | web: apexscientific.ie/gilson
Tel Dublin: 01 6854686 | Tel Cork: 021 2429776





microman[®]

Highest precision for viscous, volatile, hot or cold liquids

The new MICROMAN E positive displacement pipette offers the highest precision when pipetting problem liquids such as oil, cosmetics, blood, glycerol, radioactive samples, RNA, DNA and more!

Working like a syringe, capillary pistons eliminate the air cushion between the sample and the disposable piston, providing 100% protection against pipette contamination. With no air cushion to expand or contract, MICROMAN E retains superior accuracy and consistency when handling dense or viscous samples.

ADVANCED NEW DESIGN TO INCREASE PIPETTING PERFORMANCE



Larger Colour-Coded Push Button

More contact between thumb and button to optimise pipetting action, enhancing comfort. Same colour coding between push button and CPs



Volume Control Button

A new unique feature that secures the volume when pipetting, preventing accidental changes. Visible volume adjuster on the front of the body



New Finger Hook

Thumb is positioned higher, reducing fatigue during repetitive pipetting tasks



New Body Shape

Larger and more comfortable body to make it easier to hold in the hand



QuickSnap CP Fitting System

Unique patented feature making fitting the disposable capillary pistons as easy as fitting a regular tip

