



MICROMAN E

Ergonomic Pipetting For Problem Liquids

Guarantee accuracy, precision and comfort when pipetting volatile and viscous samples. MICROMAN E together with disposable Capillary Pistons ensures contamination-free pipetting and protects your pipette from corrosion. Get 30% off list price.

**30%
OFF**

- Experience complete and rapid pipetting of viscous, dense, hot or cold samples.
- Avoid leaks when pipetting volatile (high vapour pressure) liquids.
- No cross-contamination or residual sample carry-over.

MICROMAN E		
MODELS VOLUME	ORDER REF	OFFER PRICE
M10E 1-10 µL	FD10001	€310.00 €217.00
M25E 3-25 µL	FD10002	
M50E 20-50 µL	FD10003	
M100E 10-100 µL	FD10004	
M250E 50-250 µL	FD10005	
M1000E 100-1000 µL	FD10006	

Not to be used in conjunction with any other offer or contract discount.
Valid until 30/04/17 to ROI customers only.

How To Order:

Call Dublin: 01 6854686 or Cork: 021 2429776

Email: gilson@apexscientific.ie

Visit: www.apexscientific.ie/gilson





CAPILLARY PISTONS (CPs)

100% Protection Against Contamination

MICROMAN Capillary Pistons (CPs) offer 100% protection against contamination. Available as standard CPs or sterile CPs which are free of DNA, RNA, DNAses, RNAses, ATP, pyrogens and trace metals. Now with 30% off list price and new extended Sterile CP range.

**30%
OFF CPs**

- Large volume range 1µL to 1000 µL across six different models.
- Free of dyes, no risk of contamination from CP's components.
- Lightweight and slim CPs that fit in even the narrowest and deepest tubes.

CAPILLARY PISTONS				
CAPILLARY PISTON	ORDER REF	QTY/BOX	LIST PRICE	OFFER PRICE
CP10	F148412	192	€75.00	€52.50
	F148312	960	€296.25	€207.38
CP25	F148112	200	€62.50	€43.75
	F148012	960	€326.25	€228.38
CP50	F148113	200	€62.50	€43.75
	F148013	960	€326.25	€228.38
CP100	F148414	192	€75.00	€52.50
	F148314	960	€296.25	€207.38
CP250	F148114	200	€62.50	€43.75
	F148014	960	€326.25	€228.38
CP1000	F148560	182	€75.00	€52.50
STERILISED CAPILLARY PISTONS				
CP10ST	F148413	192	€85.00	€59.50
	F148313	960	€421.25	€294.88
CP25ST	F148712	576	€252.50	€176.75
CP50ST	F148713	576	€252.50	€176.75
CP100ST	F148415	192	€85.00	€59.50
	F148315	960	€421.25	€294.88
CP250ST	F148714	576	€252.50	€176.75
CP1000ST	F148180	182	€85.00	€59.50



Not to be used in conjunction with any other offer or contract discount.
Valid until 30/04/17 to ROI customers only.

How To Order:

Call Dublin: 01 6854686 or Cork: 021 2429776

Email: gilson@apexscientific.ie

Visit: www.apexscientific.ie/gilson